



WALK-IN INDUSTRIAL BURN-IN CHAMBERS



LABORATORY OVENS



Modular-Custom Design configurations built-to requirements



Curing, drying, ageing, annealing or any heat processing tests at higher temperatures Envisys has the oven to meet your needs


Whether it is a standard conventional drying application using shelves or conveyor drying or drying through infrared bulbs Envisys provides all types of drying solutions with custom configured designs and sizes

Parameters/Models	Specifications/Work-Space Dimensions
EO150-250T	550X400X700 MM
EO300-250T	650X600X800 MM
EO600-250T	800X850X900 MM
EO950-250T	1000X100X950 MM
EO1400-250T	1100X1350X950 MM
Volume (in liters)	150,300,600,950&1400
Temperature	250 °C



CONTROL FEATURES

Microprocessor based PID Single Set Point or Programmable controller with precise temperature control External timer for cut off of operation (Optional) Blind High temperature safety cut-off within the unit Solid state control relays for heating power control. Etc.,



ENVISYS TECHNOLOGIES PVT LTD

ISO 9001-2015 Certified Company
 #SB-118, 3rd Cross Road, I-Stage, Peenya Industrial Estate,
 Bengaluru – Karnataka – INDIA
 Land Phone : +91 80 28377888
 Fax : +91 80 28372927
 E-mail : help@envisystech.com
 Website : www.envisystech.com